




Combination Polymer Base Solid Phase Extraction

InertSep Series

- Offering high recovery rate, low blank, consuming less amount of eluent, High extraction efficiency and easy-eluting.
- Best suited for environmental water analysis, such as high volume concentration.



InertSep RP-1...MA-SDVB Acid Group Introduction Polymer.

High Flow Type, Excellent for samples high in viscosity.

InertSep MC-1...Sulfone Group Introduction Strong Cation Exchange Polymer.

InertSep MA-1...Quaternary Ammonium Group Introduction Strong Anion Exchange Polymer.

InertSep ME-1...Imino Acid Group Introduction Polymer

Selectively retains divalent cation samples.

InertSep Pharma...Nitrogen-containing-DVD Polymer.



Syringe/Barrel Types



Mini Types

•Upon Request!
Other formats/types can also be supplied.
Please contact your local distributor or
GL Sciences Inc. Japan.

InertSep Series Price List

InertSep RP-1

Description	Size	Package	Cat. No.	Price
InertSep RP-1	250mg/6mL	30EA./PKG.	5010-27000	
InertSep RP-1	30mg/1mL	100EA./PKG.	5010-27001	
InertSep RP-1	60mg/3mL	100EA./PKG.	5010-27002	
InertSep mini RP-1	230mg	50EA./PKG.	5010-27200	

InertSep MA-1

Description	Size	Package	Cat. No.	Price
InertSep MA-1	100mg/3mL	50EA./PKG.	5010-27300	
InertSep MA-1	250mg/6mL	30EA./PKG.	5010-27301	
InertSep MA-1	500mg/6mL	30EA./PKG.	5010-27302	
InertSep mini MA-1	280mg	50EA./PKG.	5010-27205	

InertSep MC-1

Description	Size	Package	Cat. No.	Price
InertSep MC-1	100mg/3mL	50EA./PKG.	5010-27350	
InertSep MC-1	250mg/6mL	30EA./PKG.	5010-27351	
InertSep MC-1	500mg/6mL	30EA./PKG.	5010-27352	
InertSep mini MC-1	280mg	50EA./PKG.	5010-27210	

InertSep Pharma

Description	Size	Package	Cat. No.	Price
InertSep Pharma	30mg/1mL	100EA./PKG.	5010-27100	
InertSep Pharma	60mg/3mL	100EA./PKG.	5010-27101	

* The specification and the cartridge types are subject to change without notice due to continual improvements.

* All brand names and product names are trademarks or registered trademarks of GL Sciences Inc.

* TM or ® are not described in this brochure.